

# LCGC<sup>®</sup>

e u r o p e

solutions for separation scientists

May 2023 | Volume 36 Number 5

[www.chromatographyonline.com](http://www.chromatographyonline.com)

## LIQUID CHROMATOGRAPHY

Improving your separation for (relatively) simple mixtures

## GAS CHROMATOGRAPHY

The impact of green chemistry on GC

## PEER REVIEW

Suspect screening of phytotoxins in the environment

## CHROMATOGRAPHY TECHNOLOGY

A review of the latest LC products

## ENVIRONMENTAL ANALYSIS FOCUS

Monitoring CECs in river water

## PRODUCTS

The latest product updates



# Source Supported Suspect Screening

A modern method to analyze phytotoxins in soil and water



# Suspect Screening of Lupin-Produced Phytotoxins in Environmental Samples

Xiaomeng Liang<sup>1</sup>, Jan H. Christensen<sup>1</sup>, Thomas D. Bucheli<sup>2</sup>, and Nikoline Juul Nielsen<sup>1</sup>, <sup>1</sup>Department of Plant and Environmental Sciences, University of Copenhagen, Frederiksberg C, Denmark, <sup>2</sup>Environmental Analytics, Agroscope, Zürich, Switzerland

Phytotoxins, the secondary metabolites synthesized by plants, are now recognized as a new category of environmental micropollutants. So far, only a limited number of phytotoxins have been detected and reported in terrestrial and aquatic environments, partially due to the analytical challenge involved. In this article, a robust, reliable, and efficient approach, namely *Source Supported Suspect Screening (4S)*, is introduced for high-throughput analysis of phytotoxins from their plant origin to the downstream environmental compartments. The approach was established on an ultrahigh-performance liquid chromatography–electrospray ionization–quadrupole-time-of-flight-mass spectrometry (UHPLC–ESI–QTOF-MS) analytical platform where optimal sample preparation and chromatographic conditions were investigated and a new suspect screening protocol was suggested.

## KEY POINTS

- Selective pressurized liquid extraction and multilayer solid-phase extraction methods are described for high-throughput plant, soil, and water sample preparations.
- Optimal analytical conditions are described for the improvement of detection sensitivities and coverage.
- A Source Supported Suspect Screening (4S) approach is described; phytotoxins detected in the source plant were used to improve the identification of phytotoxins in soil and water.

Phytotoxins are secondary metabolites, typically <1000 Da, synthesized by plants. The study of phytotoxins as endogenous plant hormones, dietary supplements, cosmetic ingredients, and pharmaceuticals dates back decades and is still of interest today. Most recently, evidence of phytotoxins transferring from plant to terrestrial and aquatic environments has been reported, which demonstrates that these natural toxins are a new type of micropollutant in the environment (1–4). One example is Balkan endemic nephropathy, a disease that mostly occurs in the rural Balkan regions, and is linked to aristolochic acids, the nephrotoxic phytotoxins produced from *Aristolochia clematitis* L.. In this case, aristolochic acids were reported as persistent groundwater contaminants in Serbia, with groundwater being the major source of cooking and drinking water for the local residents (1).

Analyzing phytotoxins in environmental matrices such as soil, groundwater, and surface water is challenging. The compounds are structurally highly diverse. Their environmental occurrences and accumulations can be influenced by many factors such as vegetation, geographical conditions, agronomic factors, as well as season, weather, and climate change (5,6). In addition, the lack of reference standards and informative databases makes it difficult to perform qualitative analysis (screening and compound identification) and quantitative analysis (concentration measurement).

A robust, reliable, and efficient approach, namely *Source Supported Suspect Screening (4S)*, has been introduced for high-throughput analysis of a wide range of phytotoxins, from their plant origin to the downstream environmental compartments, such as soil, groundwater, and surface water. The approach was demonstrated in a field study of *Lupinus angustifolius* L. (also known as *blue lupin*). This plant is a promising alternative to soya beans in European farming systems to provide food and feed; it is also a popular green manure in marginal lands (7). A limited number of lupin-derived phototoxic compounds, mainly alkaloids and flavonoids,

have been reported in agricultural soil and soil pore water (8,9). In this study, it was hypothesized that there could be a much higher number of phytotoxins transferred from the lupin plant to the surrounding soil and water environments. Based on this hypothesis, 106 plant, soil, drainage, and stream water samples were collected on a regular basis before, during, and after the growth cycle. These were analyzed through a 5-month agriculture cultivation of blue lupin in northern Switzerland. This work has been reported previously in the journal of *Environmental Science & Technology* (10).

### Experimental and Discussion Sample Preparation for Plant, Soil, and Water:

In the study, 41 lupin plant samples and 26 root soil samples were collected and extracted following a previously reported selective pressurized liquid extraction (sPLE) method (11). Extractions were performed on an ASE (accelerated solvent extraction) instrument (Dionex Corporation) under high pressure (1500 psi). The high pressure allowed the liquid extraction (50% aqueous methanol) to be performed at an elevated temperature (100 °C) with high recovery yields (67–103% for eight reference analytes, that is, convallatoxin,  $\kappa$ -strophanthidin, digoxin, digitoxigenin, odoroside A, proscillaridin, bufalin, and withanolide A) (11). In comparison to other extraction methods, such as Soxhlet and ultrasonication extractions, sPLE has several advantages, including being programmable, automated, protective to light-sensitive and easily oxidized analytes, and in-cell cleanup. The sorbent layer (C18) on the bottom of the extraction cell not only adds selectivity but also enables the extraction-filtration cleanup to be achieved in one step. Key findings in the sPLE study (11) showed that the conventional heating step prior to solvent introduction

caused significant thermal degradation (observed as reduced recovery) of a specific phytotoxin, digoxin, whereas a heating step after solvent introduction did not. Furthermore, matrix effects in the subsequent liquid chromatography–electrospray ionization–mass spectrometry (LC–ESI–MS) analysis were highly dependent on the sPLE dispersant used

(Ottawa sand, diatomaceous earth, or glass beads), with the use of glass beads resulting in insignificant matrix effects. The matrix effects were also dependent on the combination of solvent strength and extraction temperature; the selection of 50% aqueous methanol and 100 °C kept matrix effects to a minimum, with an absolute average of 10% ( $\sigma = 4\%$ ).

**YMC**  
EUROPE GMBH  
The Selectivity Company

**HPLC 2023**  
DÜSSELDORF  
18-22 June 2023

Visit us at  
booth #6!


YMC HPLC COLUMN  
notes:  
YMC Co., Ltd. Japan

**Experts in  
Reproducibility**


- **Robust Bio-RP (U)HPLC**  
Extremely inert particles produce sharp peaks for proteins/peptides, oligonucleotides and mAbs.
- **High Recovery IEX**  
Low adsorption and excellent resolution in analysis of proteins, mAbs and oligonucleotides.
- **Highly Efficient HIC & SEC**  
Different selectivities for fast and reliable analysis of proteins, mAbs and ADCs.



**FIGURE 1:** Extracted ion chromatograms of digitoxigenin, bufalin, genistein, and artemisinin. Column temperature was set at 25, 35, and 45 °C. The mobile phase was (A) water and (B) 95% methanol, both with 5 mM ammonium formate as additive and prepared at pH 3.0. Flow rate was 0.3 mL/min. The UHPLC program was 0–1 min 1% B followed by linear gradient to 100% B in 10 min. The y-axes are linked for each of the phytotoxins. Adapted and reproduced with permission from reference 15.



**FIGURE 2:** Flow-diagram of the Source Supported Suspect Screening (4S) approach. Phytotoxins detected in the source plant were used to improve the identification of phytotoxins in soil and water. Adapted and reproduced with permission from reference 10.



In parallel, 30 drainage water and nine stream water samples were collected from the field itself (drainage water) and downstream from the field (stream water). These samples were concentrated following an established multilayer solid-phase extraction (mSPE) method reported by Mechelke et al. (12). The four-layer sorbents provided a broad coverage of analytes ranging from highly polar to nonpolar ( $\log D_{\text{ow,pH7}}^{-14-8}$ ), with overall good recovery (> 70%) and low matrix effects (63–72% on average), according to Mechelke et al. and Günthardt et al. (12,13).

**Instrumental Analysis:** The analysis was performed on an ultrahigh-performance liquid chromatography–electrospray ionization–quadrupole–time-of-flight–mass spectrometry<sup>E</sup> (UHPLC–ESI–QTOF–MS<sup>E</sup>) analytical platform (Waters). A reverse-phase-based charged surface hybrid (CSH) column was chosen for its wide selectivity and good peak shape for basic compounds (14). Prior to batch analysis, chromatographic factors such as mobile phase pH/additive, organic modifier, column temperature, and gradient time were evaluated (15). Results showed 5 mM ammonium formate

(pH 6.5) as mobile phase additive and methanol as organic modifier significantly improved the overall detection sensitivity compared to the popular 0.1% formic acid and acetonitrile recipe. The optimal column temperature was set to 25 °C to avoid thermal degradation of phytotoxins during the chromatographic separation (Figure 1). The optimal gradient time was 30 min as a balance of peak capacity, peak production rate, and the total analysis time. ESI<sup>+</sup> has been shown to efficiently ionize alkaloids, flavonoids, and many other phytotoxins and was therefore selected as the ionization mode (15). The LC–MS data were acquired using the data-independent-acquisition (DIA) mode where all precursor ions were fragmented without pre-selection. DIA is advantageous in terms of higher sample coverage and better reproducibility between samples, allowing the detection, structural characterization, and quantification to be achieved as an “all-in-one” analytical run (16).

### Suspect Screening by the 4S Approach

For suspect screening, a list containing 733 phytotoxins was created. The list was a combination of 516 prioritized phytotoxins from the Toxic Plants–PhytoToxins database (17) and 217 phytotoxins from literature reports of secondary metabolites derived from *Lupinus* spp. (18–20). As demonstrated in Figure 2, raw LC–MS data were divided into two datasets: (i) a source plant dataset and (ii) an environmental fates (soil/water) dataset. A suspect screening was first performed on the plant dataset, which yielded a list of 66 tentatively identified phytotoxins in the lupin plant. The initial suspect screening list was subsequently updated by including the LC retention time of phytotoxins identified with high confidence (level 2) and phytotoxins confirmed by reference standards



(level 1). The second suspect screening was performed on the soil/water dataset, where the updated suspect screening list was used. The retention times obtained from the plant dataset were included together with MS/MS fragmentation patterns as one of the level 2 identification criteria for phytotoxin candidates occurring in the soil/water matrix.

### Phytotoxins in Plant, Soil, and Water

Using the 4S approach, 72 phytotoxins from the lupin plant and the environmental samples were identified with high confidence (level 1 and 2). As shown in Figure 2, they were exclusively flavonoids and alkaloids, mainly genistein, apigenin, chrysoeriol, wighteone, luteone, lupanine, sparteine, angustifoline, gramine, ammodendrine, formononetin, daidzein, caffeine, senecionine, retrorsine N-oxide, chelidonine, and/or their derivatives. The identities of genistin, genistein, gramine, sparteine, formononetin, daidzein, caffeine, senecionine, retrorsine N-oxide, and chelidonine were confirmed by comparison to reference standards. Sixty-six phytotoxins were consistently detected in the plant over its growing season. Of the lupin-derived phytotoxins, 45 were detected in the soil, 34 were detected in drainage water, and 13 were detected in stream water. The results demonstrated that phytotoxins may be released from lupin and transferred to the surrounding terrestrial and aquatic environment. The established 4S approach can be powerful to assist in the discovery of phytotoxins in the environment.

### Conclusion

As revealed in the study, the established 4S approach can be an efficient and reliable strategy to detect phytotoxins in multiple matrices, and can backtrack their occurrences from the downstream

environments to the plant origin. In the long term, we expect the approach to assist in the discovery of more terrestrial- and aquatic-occurring phytotoxins, remedying the knowledge gap between their origin and environmental fate, and initiating related risk assessments and safety regulation.

### References

- 1) Tung, K.-K.; Chan, C.-K.; Zhao, Y.; et al. Occurrence and Environmental Stability of Aristolochic Acids in Groundwater Collected from Serbia: Links to Human Exposure and Balkan Endemic Nephropathy. *Environ. Sci. Technol.* **2019**, *54* (3), 1554–1561. DOI: DOI:10.1021/acs.est.9b05337



### SEC-MALS/MS FOR ACCURATE SAMPLE CHARACTERIZATION



#### A column optimized for MALS/MS and biopharmaceutical analysis

The **TSKgel® UP-SW3000-LS** U/HPLC SEC column offers unique noise suppression, resulting in increased sensitivity of advanced detection.



#### A MALS detector with revolutionary technology

The **LenS<sub>3</sub>®** MALS detector's design eliminates noise from stray light, thereby maximizing S/N. This results in incredibly sensitive and accurate biomolecular MW measurements.



#### A team of experts to support your work

Our team of chromatography experts provides our biopharma partners with solutions to develop safe and efficient therapies.



Tosoh Bioscience and TSKgel are registered trademarks of Tosoh Corporation. LenS is a registered trademark of Tosoh Bioscience LLC in the USA, India, and Japan.

Join us at our Booth G1 at HPLC, June 18-22 in Düsseldorf, Germany!

+ 49 6155 7043700

info.tbg@tosoh.com

www.tosohbioscience.de

**TOSOH BIOSCIENCE**

- 2) Picardo, M.; Sanchís, J.; Núñez, O.; Farré, M. Suspect Screening of Natural Toxins in Surface and Drinking Water by High Performance Liquid Chromatography and High-Resolution Mass Spectrometry. *Chemosphere* **2020**, *261*, 127888. DOI: 10.1016/j.chemosphere.2020.127888
- 3) Nanusha, M. Y.; Krauss, M.; Brack, W. Non-Target Screening for Detecting the Occurrence of Plant Metabolites in River Waters. *Environ. Sci. Eur.* **2020**, *32* (1), 1–14. DOI: 10.1186/s12302-020-00415-5
- 4) Hansen, H. C. B.; Hilscherova, K.; Bucheli, T. D. Natural Toxins: Environmental Contaminants Calling for Attention. *Environ. Sci. Eur.* **2021**, *33* (1), 112. DOI: 10.1186/s12302-021-00543-6
- 5) Bucheli, T. D. Phytotoxins: Environmental Micropollutants of Concern? *Environ. Sci. Tech.* **2014**, *48* (22), 13027–13033. DOI: 10.1021/es504342w
- 6) Günthardt, B. F.; Hollender, J.; Scheringer, M.; et al. Aquatic Occurrence of Phytotoxins in Small Streams Triggered by Biogeography, Vegetation Growth Stage, and Precipitation. *Sci. Total Environ.* **2021**, *798*, 149128. DOI: 10.1016/j.scitotenv.2021.149128
- 7) Lopez-Bellido, L.; Fuente, M. Lupin Crop as an Alternative Source of Protein. *Adv. Agron.* **1986**, *40*, 239–295. DOI: 10.1016/S0065-2113(08)60284-9
- 8) Hama, J. R.; Jorgensen, D. B. G.; Diamantopoulos, E.; et al. Indole and Quinolizidine Alkaloids from Blue Lupin Leach to Agricultural Drainage Water. *Sci. Total Environ.* **2022**, *834*, 155283. DOI: 10.1016/j.scitotenv.2022.155283
- 9) Andersen, I. K.; Laursen, B. B.; Rasmussen, J.; Fomsgaard, I. S. Optimised Extraction and LC-MS/MS Analysis of Flavonoids Reveal Large Field Variation in Exudation into *Lupinus Angustifolius* L. Rhizosphere Soil. *Rhizosphere* **2022**, *22* (4), 100516. DOI:10.1016/j.rhisph.2022.100516
- 10) Liang, X.; Christensen, J. H.; Bucheli, T. D.; Nielsen, N. J. Source-Supported Suspect Screening (4S) of Phytotoxins in Terrestrial and Aquatic Environments: A Field Study of *Lupinus Angustifolius* L. (Blue Lupin). *Environ. Sci. Technol.* **2023**, *57* (6), 2333–2340. DOI: 10.1021/acs.est.2c05387
- 11) Liang, X.; Nielsen, N. J.; Christensen, J. H. Selective Pressurized Liquid Extraction of Plant Secondary Metabolites: *Convallaria Majalis* L. as a Case. *Anal. Chim. Acta: X* **2020**, *4*, 100040. DOI: 10.1016/j.acax.2020.100040
- 12) Mechelke, J.; Longrée, P.; Singer, H.; Hollender, J. Vacuum-Assisted Evaporative Concentration Combined with LC-HRMS/MS for Ultra-Trace-Level Screening of Organic Micropollutants in Environmental Water Samples. *Anal. Bioanal. Chem.* **2019**, *411*, 2555–2567. DOI: 10.1007/s00216-019-01696-3
- 13) Günthardt, B. F.; Schönsee, C. D.; Hollender, J.; et al. "Is There Anybody Else Out There?"—First Insights from a Suspect Screening for Phytotoxins in Surface Water. *Chimia* **2020**, *74* (3), 129–135. DOI: 10.1016/j.scitotenv.2021.149128
- 14) Fountain, K. J.; Hewitson, H. B.; Iraneta, P. C.; Morrison, D. Practical Applications of Charged Surface Hybrid (CSH) Technology. Waters Application Note 720003720EN, 2010.
- 15) Liang, X.; Christensen, J. H.; Nielsen, N. J. Enhancing the Power of Liquid Chromatography–Mass Spectrometry for Chemical Fingerprinting of Phytotoxins in the Environment. *J. Chromatogr. A* **2021**, *1642*, 462027. DOI: 10.1016/j.chroma.2021.462027
- 16) Wang, R.; Yin, Y.; Zhu, Z. -J. Advancing Untargeted Metabolomics Using Data-Independent Acquisition Mass Spectrometry Technology. *Anal. Bioanal. Chem.* **2019**, *411*, 4349–4357. DOI: 10.1007/s00216-019-01709-1
- 17) Günthardt, B. F.; Hollender, J.; Hungerbühler, K.; Scheringer, M.; Bucheli, T. D. Comprehensive Toxic Plants–Phytotoxins Database and Its Application in Assessing Aquatic Micropollution Potential. *J. Agr. Food Chem.* **2018**, *66* (29), 7577–7588. DOI: 10.1021/acs.jafc.8b01639
- 18) Wink, M.; Meißner, C.; Witte, L. Patterns of Quinolizidine Alkaloids in 56 Species of the Genus *Lupinus*. *Phytochemistry* **1995**, *38* (1), 139–153. DOI: 10.1016/0031-9422(95)91890-D
- 19) Erdemoglu, N.; Ozkan, S.; Tosun, F. Alkaloid Profile and Antimicrobial Activity of *Lupinus Angustifolius* L. Alkaloid Extract. *Phytochemistry Rev.* **2007**, *6* (1), 197–201. DOI: 10.1007/s11101-006-9055-8
- 20) Von Baer, D.; Saelzer, R.; Vega, M.; et al. Isoflavones in *Lupinus Albus* and *Lupinus Angustifolius*: Quantitative Determination by Capillary Zone Electrophoresis, Evolution of Their Concentration During Plant Development and Effect on Anthracnose Causing Fungus *Colletotrichum Lupini*. *J. Chil. Chem. Soc.* **2006**, *51* (4), 1025–1029. DOI: 10.4067/S0717-97072006000400006

## CHIRALPAK® IM NOW AVAILABLE IN 3 µm

Backed by 40+ years of experience and innovation, Daicel Chiral Technologies continues to invest in your project needs

CHIRALPAK® IM is the newest in our line of robust, versatile immobilised phases

- Robust to withstand all mobile phase combinations
- Designed for HPLC and SFC to improve methods for challenging separations
- Available in 3 and 5 µm particle size

The best separations demand the best columns.

Visit Us at HPLC 2023 - Foyer #5  
Duesseldorf, Germany  
June 18 - 22, 2023

Visit us at [chiraltech.com](http://chiraltech.com)

**DAICEL**  
Chiral Technologies

**Xiaomeng Liang** was a PhD student in the Department of Plant and Environmental Sciences at the University of Copenhagen, Denmark.  
**Jan H. Christensen** is a professor in the Department of Plant and Environmental Sciences at the University of Copenhagen, and a member of the EAB for LCGC Europe.  
**Thomas D. Bucheli** leads the research group Environmental Analytics at Agroscope in Switzerland.  
**Nikoline Juul Nielsen** is an associate professor in the Department of Plant and Environmental Sciences at the University of Copenhagen.